

**Section 1. Product and Company Identification**

Item Number.: s2120-1  
 Common Name.: Oil Red O Solution  
 Intended Use : In Vitro Diagnostic use. Laboratory Use Only  
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer.: Poly Scientific R&D Corp.

70 Cleveland Ave  
 Bay Shore NY 11706

polyrnd@polyrnd.com

**Section 2. Hazard Identification**

225 Flammable Liquids Cat 2  
 320 Serious eye damage/eye irritation Cat 2B  
 336 Specific target organ toxicity, single exposure; Narcotic effects Cat 3



Danger

Highly flammable liquid and vapour.  
 Causes eye irritation.

May cause drowsiness and dizziness.

Keep away from heat/sparks/open flames/hot surfaces. No smoking. Keep container tightly closed. Ground/bond container and receiving equipment. Use explosion-proof electrical/ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Avoid breathing dust/fume/gas/mist/vapours/spray. Wash hands/skin thoroughly after handling. Use only in well-ventilated area. Wear protective gloves/protective clothing/eye protection/face protection. . In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Rinse mouth. Do not induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF eye irritation persists: Get medical advice/attention. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician. Store locked up in a well ventilated place and keep cool. Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

**Section 3. Composition Information**

Exposure Limits(A blank value indicates no information available) The specific chemical identity and/or exact percentage (concentration) of composition has been withheld as a trade secret.

Component	CAS#	PEL(mg/m3)	STEL(mg/m3)	CEIL(mg/m3)	Concentration Range
Ethyl Alcohol	64-17-5		1,900.00		25-50%
Isopropyl Alcohol	67-63-0	1,225.00	980.00		<5%
Acetone	67-64-1	1,000.00	500.00		25-50%
Methyl Alcohol	67-56-1	325.00	260.00		<5%

**Section 4. First Aid Measures**

Eye Contact : Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.

Skin Contact : Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention

Ingestion : Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention

Inhalation : Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention

The most important known symptoms and effects are described in section 2 and/or section 11.

**Section 5. Fire Fighting Measures**

Extinguishing Media.: Use Dry Chemical, Foam or Carbon Dioxide

Special Fire and Explosion Remarks : N/A

**Section 6. Accidental Release Measures**

Spill Cleanup and Disposal Special : Eliminate all ignition sources. Pick up with absorbent material.

Spill Cleanup.: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

**Section 7. Handling and Storage**

Storage and Handling Special.: N/A

Storage and handling.: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

**Section 8. Exposure Controls/Personal Protection**

Personal Protective Equipment : safety Glasses, Gloves, Vapor Respirator

This information is provided as a guide but proper PPE can only be determined by the end user and their situation.

Engineering Controls : Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

**Section 9. Physical and Chemical Properties**

Appearance.....: Turbid red liquid	Evaporation Rate.....: N/A	Water Soluable?.....: Yes
Odor.....: Sharp, penetrating	Upper Flammability Limit (%): N/A	Volatile Percent.....: 100
Odor Threshold.....: N/A	Lower Flammability Limit (%): N/A	Partition Coefficient.....: n-octanol/water: N/A
pH.....: N/A	Specific Gravity (@20C).....: 07915	Auto Ignition Temp.....: N/A
Melting Point.....: N/A	Vapor Pressure (mm Hg).....: 184	Decomposition Temp.....: N/A
Boiling Point.....: N/A	Vapor Density (Air=1).....: 2	Viscosity.....: N/A
Flash Point (F) TCC.....: N/A	Relative Density.....: N/A	

**Section 10. Stability and reactivity**

Special Remarks on Stability.: Stable

Special Remarks on Reactivity.: N/A

Water Reactive.: No

**Section 11. Toxicological Information**

Routes of Entry.: Inhalation, Skin Absorption, Ingestion

Animal Toxicity.: Ethyl Alcohol: Acute Oral (LD50) 7060 mg/kg (Rat) Acute Dermal (LD50) 500mg/24hr (Severe) (Rabbit); Methyl Alcohol: Acute Oral (LD50) 5628 mg/kg (Rat) Acute Dermal (LD50) 500mg/24hr (Mod Irrit) (Rabbit); Isopropyl Alcohol: Acute Oral (LD50) 5045 mg/kg (Rat); Oil Red O: N/A; Acetone: Acute Oral (LD50) 9.75 gm/kg (Rat) Acute Dermal (LD50) 20

Human Toxic Effects.: Target Organs: Respiratory system, skin, eyes, CNS, liver, blood, Kidney, Reprod System, Liver, CNS

Potential Acute Health Effects...: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects...: Isopropyl Alcohol: IARC Code 3

**Section 12. Ecological Information**

Ecological Information.: N/A

**Section 13. Disposal Considerations**

Waste Disposal.: Dispose of in accordance with local, state and federal laws.

**Section 14. Transport Information**

DOT Identification.: UN1090, Acetone,3,II

**Section 15. Regulatory Information**

State Regulations.: New York release reporting list: Acetone

Component	CAS#	Sara Section 311 Reporting									
		Acute	Chronic	Fire	Pressure	Reactive	SARA302	SARA313	CERCLA	RCRA	
Ethyl Alcohol	64-17-5	No	No	No	No	No	Yes	Yes	No		
Isopropyl Alcohol	67-63-0	No	No	No	No	No	Yes	No	No		
Acetone	67-64-1	No	No	No	No	No	No	Yes	No		
Methyl Alcohol	67-56-1	No	No	No	No	No	Yes	Yes	No		

**Section 16. Other Information**

Review Date : 3/15/2023

Reviewed by : Admin

MSDS Group Id.: 190

Notice: This SDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this SDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.