

Section 1. Product and Company Identification

Item Number.: s260-1
 Common Name.: Poly Mount (Toluene Based)
 Intended Use : In Vitro Diagnostic use. Laboratory Use Only
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer.: Poly Scientific R&D Corp.

70 Cleveland Ave
 Bay Shore NY 11706

polyrnd@polyrnd.com

Section 2. Hazard Identification

225 Flammable Liquids Cat 2
 304 Aspiration hazard Cat 1
 315 Skin corrosion/irritation Cat 2
 336 Specific target organ toxicity, single exposure; Narcotic effects Cat 3
 361 Reproductive toxicity Cat 2
 370 Specific target organ toxicity, single exposure Cat 1

**Danger**

Highly flammable liquid and vapor. May be fatal if swallowed and enters airways. Causes skin irritation. May cause drowsiness or dizziness. Suspected of damaging the unborn child. Causes damage to CNS, Kidneys and Liver.

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Keep away from heat/sparks/open flames/hot surfaces. - No smoking. Use explosion-proof electrical/ventilating/lighting equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Do not breathe dust/fume/gas/mist/vapours/spray. Wash hands/skin thoroughly after handling. Use in a well-ventilated area. Use personal protective equipment as required. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN (or hair): Take off contaminated clothing. Wash skin with soapy water/shower. If skin irritation or rash occurs: Get medical advice/attention. Wash contaminated clothing before reuse. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Call a poison center/doctor if you feel unwell. Immediately call a POISON CENTER or doctor/physician. If exposed or concerned: Get medical advice/attention. In case of fire: Use appropriate media to extinguish. Store locked up. Disposal Dispose of contents/container in accordance with local/regional/national/international regulations. Store in a cool, well-ventilated place, tightly closed and locked up.

Hazards not otherwise covered by GHS: None

Section 3. Composition Information

Exposure Limits(A blank value indicates no information available) The specific chemical identity and/or exact percentage (concentration) of composition has been withheld as a trade secret.

Component	CAS#	PEL(mg/m3)	STEL(mg/m3)	CEIL(mg/m3)	Concentration Range
Tricresyl Phosphate	1330-78-5		0.10		<5%
Toluene	108-88-3		188.00		25-50%

Section 4. First Aid Measures

Eye Contact : Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.
 Skin Contact : Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention
 Ingestion : Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention
 Inhalation : Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention
 The most important known symptoms and effects are described in section 2 and/or section 11.

Section 5. Fire Fighting Measures

Extinguishing Media.: Dry Chemical, Carbon Dioxide, Water Spray, Alcohol Foam
 Special Fire and Explosion Remarks...: N/A

Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special...: Pick up spill with absorbent vermiculite
 Spill Cleanup.: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

Section 7. Handling and Storage

Storage and Handling Special...: N/A
 Storage and handling.: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment...: Splash Goggles, Gloves, Synthetic Apron, vapor Respirator
 This information is provided as a guide but proper PPE can only be determined by the end user and their situation.
 Engineering Controls.: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

Section 9. Physical and Chemical Properties

Appearance.....: Colorless liquid	Evaporation Rate.....: N/A	Water Soluable?.....: Yes
Odor.....: Benzene like	Upper Flammability Limit (%): N/A	Volatile Percent.....: 100
Odor Threshold.....: N/A	Lower Flammability Limit (%): N/A	Partition Coefficient.....: n-octanol/water: N/A
pH.....: NA	Specific Gravity (@20C).....: 0.87	Auto Ignition Temp.....: N/A
Melting Point.....: N/A	Vapor Pressure (mm Hg).....: 53	Decomposition Temp.....: N/A
Boiling Point.....: N/A	Vapor Density (Air=1).....: 3.17	Viscosity.....: N/A
Flash Point (F) TCC.....: N/A	Relative Density.....: N/A	

Section 10. Stability and reactivity

Special Remarks on Stability...: Stable

Special Remarks on Reactivity...: NA

Water Reactive.: No

Section 11. Toxicological Information

Routes of Entry.: Inhalation, Skin Absorption, Ingestion

Animal Toxicity.: Toluene: Acute Oral (LD50) 5000 mg/kg (Rat) I.V. (LD50) 1960 mg/kg (Rat); Elvacite: N/A; Tricresyl Phosphate: Acute Oral (LD50) 5190 mg/kg (Rat) Acute Dermal (LD50) 1500 mg/kg (Cat)

Human Toxic Effects.: Target Organs: Kidneys, Brain, Bladder, Liver

Potential Acute Health Effects...: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects...: Toluene: IARC Code 3

Section 12. Ecological Information

Ecological Information.: N/A

Section 13. Disposal Considerations

Waste Disposal.: Dispose of in accordance with local, state and federal laws.

Section 14. Transport Information

DOT Identification.: UN1294, Toluene, 3,II

Section 15. Regulatory Information

State Regulations.: New York release reporting list: Toluene

Component	CAS#	Sara Section 311 Reporting									
		Acute	Chronic	Fire	Pressure	Reactive	SARA302	SARA313	CERCLA	RCRA	
Tricresyl Phosphate	1330-78-5	No	No	No	No	No	No	No	No	No	
Toluene	108-88-3	No	No	No	No	No	Yes	Yes	No		

Section 16. Other Information

Review Date : 3/15/2023

Reviewed by : Admin

MSDS Group Id.: 99

Notice: This SDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this SDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.